


Becher, dickwandig
(Filtrierbecher) mit
Teilung, mit Ausguss

Beaker, heavy wall
(Jars, filtering) with
graduation, with spout




**Béchers à parois
épaisses**
(béchers à filtration)
avec graduation et bec

Vasos, pared gruesa
(para filtración)
graduados, con pico

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Originalpkg./Stück Stand. pack/quantity
21 131 24	100	52	85	10
21 131 29	150	54	93	10
21 131 36	250	70	94	10
21 131 44	500	89	124	10
21 131 54	1 000	105	160	10
21 131 63	2 000	135	195	10
21 131 68	3 000	157	205	4
21 131 73	5 000	182	256	1
21 131 86 ¹	10 000	225	340	1
21 131 88 ¹	15 000	260	390	1
21 131 91 ¹	20 000	285	430	1

¹ Ohne Teilung
¹ Without graduation
¹ Sans graduation
¹ Sin graduación




Philipsbecher
mit Ausguss

Beakers, Philips
with spout

Béchers de Philips
avec bec

Vasos Philips
con pico

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Originalpkg./Stück Stand. pack/quantity
21 141 29	150	59	87	10
21 141 36	250	68	105	10
21 141 44	500	86	142	10



Berzeliusbecher
Beaker, Berzelius

Béchers de Berzelius
Vasos Berzelius

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Originalpkg./Stück Stand. pack/quantity
21 126 01	100	50	78	10