

Ansetzflaschen
(Erlenmeyerform)
weithalsig

Flacons coniques
(forme Erlenmeyer)
à col large


Conical flasks
wide neck

Frascos para disoluciones
(forma Erlenmeyer)
cuello ancho

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d ₁ mm	h mm	Originalpkg./Stück Stand. pack/quantity
21 227 68 ^{1, 2}	3 000	190	106	285	1
21 227 73 ^{1, 2}	5 000	220	108	322	1
21 227 86 ^{1, 2}	10 000	285	147	420	1

¹ Nicht nach DIN
¹ Non-DIN size
¹ Non conforme aux dimensions DIN
¹ No según DIN

² ohne Teilung
² without graduation
² sans graduation
² sin graduados



Erlenmeyer-Kolben
mit DIN-Gewinde und
Teilung

Fioles Erlenmeyer
avec filetage DIN,
et graduation



Screwthread
Erlenmeyer flasks
with DIN threads and
graduation

Matraces Erlenmeyer
con rosca DIN,
graduados

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN thread GL	d mm	h mm	Originalpkg./Stück Stand. pack/quantity
21 803 24	100	GL 25	64	105	10
21 803 24 5 ¹	100	GL 25	64	109	10
21 803 36	250	GL 32	85	145	10
21 803 36 5 ¹	250	GL 32	85	149	10
21 803 44	500	GL 32	105	175	10
21 803 44 5 ¹	500	GL 32	105	180	10
21 803 54	1 000	GL 32	131	220	10
21 803 54 5 ¹	1 000	GL 32	131	225	10

¹ Mit Schraubverschluss-Kappe, Ersatzkappen siehe Seite 128, 129.
¹ With screw-cap, closed, for spare caps please see page 128, 129.
¹ Avec capuchon vissé, capuchons vissables voir page 128, 129.
¹ Con tapa de plástico, caperuzas de unión roscada ver pág. 128, 129.

Achtung: Nicht mit fest verschlossener Kappe autoklavieren.
Important: Do not autoclave bottles with a tightly screwed cap.
Attention: ne pas autoclaver avec le bouchon vissé à fond.
Atención: no someter al autoclave con la tapa totalmente cerrada.



Jodzalkolben
mit Normschliff und
Deckelstopfen, hohl

Fioles pour détermination
de l'indice d'iode
à RIN et bouchons six
pans creux

Flasks with standard
ground joint
and with hollow, flat-
head stopper (iodine
determination)

Matraces para índice
de yodo, con esmerilado
normalizado y tapón
hueco

Best.-Nr. Cat. No.	Inhalt Capacity ml	NS Joint	d mm	h mm	Originalpkg./Stück Stand. pack/quantity
24 192 27	100	29/32	64	120	10
24 192 37	250	29/32	85	160	10
24 192 46	500	29/32	105	195	10
24 192 56	1 000	29/32	131	235	10

Ersatzstopfen: 24 140 09
Spare stopper: 24 140 09
Bouchons six pans creux: 24 140 09
Tapón hueco: 24 140 09