

Kristallisierschalen

Cristallisoirs

Crystallizing dishes

Cristalizadores


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Originalpckg./Stück Stand. pack/quantity
mit Ausguss DIN 12 338/with spout DIN 12 338/avec bec DIN 12 338/con pico DIN 12 338				
21 311 24 <sup>1</sup>	20	40	25	10
21 311 32 <sup>1</sup>	40	50	30	10
21 311 34 <sup>1</sup>	60	60	35	10
21 311 38	100	70	40	10
21 311 41	150	80	45	10
21 311 44	300	95	55	10
21 311 49	500	115	65	10
21 311 54	900	140	75	10
21 311 59	2 000	190	90	10
21 311 63	3 500	230	100	10
ohne Ausguss DIN 12337/without spout DIN 12337/sans bec DIN 12337/sin pico DIN 12337				
21 313 24 <sup>1</sup>	20	40	25	10
21 313 32 <sup>1</sup>	40	50	30	10
21 313 34 <sup>1</sup>	60	60	35	10
21 313 38	100	70	40	10
21 313 41	150	80	45	10
21 313 44	300	95	55	10
21 313 49	500	115	65	10
21 313 54	900	140	75	10
21 313 59	2 000	190	90	10
21 313 63	3 500	230	100	10

<sup>1</sup> Ohne Aufdruck

<sup>1</sup> Without print

<sup>1</sup> Sans impression

<sup>1</sup> Sin marca




Wägegläser  
mit Deckel

Pèse-filtres  
avec couvercle à bouton

Weighing bottles  
with grip stopper

Frascos para pesar  
con tapa



Best.-Nr. Cat. No.	Inhalt Capacity ml	d x h mm	Originalpckg./Stück Stand. pack/quantity
niedrige Form/flat shape/forme basse/forma baja			
24 210 13	5	28 x 25	10
24 210 23	15	38 x 30	10
24 210 32	30	54 x 30	10
24 210 41	80	85 x 30	10
hohe Form/tall shape/forme haute/forma alta			
24 211 13	10	28 x 40	10
24 211 18	20	32 x 50	10
24 211 23	45	38 x 70	10
24 211 24	70	44 x 80	10