


Best.-Nr. Cat. No.	d mm	Originalpkg./Stück Stand. pack/quantity
DURAN		
21 321 24	40	10
21 321 32	50	10
21 321 34	60	10
21 321 41	80	10
21 321 46	100	10
21 321 52	125	10
21 321 57	150	10
21 321 61	200	10
21 321 66	250	1
aus Kalk-Soda-Glas/in soda-lime glass/en verre sodo-calcique/en vidrio cal-soda		
23 321 24	40	10
23 321 32	50	10
23 321 34	60	10
23 321 38	70	10
23 321 41	80	10
23 321 43	90	10
23 321 46	100	10
23 321 51	120	10
23 321 52	125	10
23 321 57	150	10
23 321 61	200	10
23 321 66	250	1



Uhrglasschalen
Rand verschmolzen,
DIN 12 341


Verres de montre
bord rébrûlé,
DIN 12 341

Watch glass dishes
with fused edges,
DIN 12 341

Vidrios de reloj
borde quemado,
DIN 12 341

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Originalpkg./Stück Stand. pack/quantity
24 204 23	75	50	70	10
24 204 24	100	54	75	10
24 204 36 ¹	250	74	102	10

¹ Nur solange Vorrat reicht
¹ Discontinued, available until stocks exhausted
¹ Uniquement jusqu'à épuisement du stock
¹ Sólo hasta agotar existencias




Organgläser
ohne Stopfen

Verres à organes
sans bouchon

Organ storage jars
without stoppers

Frascos para órganos
sin tapón

Best.-Nr. Cat. No.	d x h mm	Volumen ca. Volume approx. ml	Originalpkg./Stück Stand. pack/quantity
24 208 34	60 x 40	75	10
24 208 41	80 x 50	175	10
24 208 45	100 x 60	325	10
24 208 57	150 x 80	1 000	10



Dosen
mit Deckel

Boîtes
avec couvercle

Jars
with lid

Tarros
con tapa