


Best.-Nr. Cat. No.	d x h mm	Volumen ca. Volume approx. ml	Originalpckg./Stück. Stand. pack/quantity
24 207 34	60 x 35	70	10
24 207 45	103 x 55	250	10
24 207 51	121 x 64	500	10

**Dosen**  
mit Falz und Deckel,  
nach DIN 12 340

**Boîtes**  
avec épaulement et  
couvrecl, DIN 12 340

**Jars**  
with shoulder and lid,  
DIN 12 340

**Tarros**  
con borde y tapa, según  
DIN 12 340



Best.-Nr. Cat. No.	d x h mm	Volumen ca. Volume approx. ml	Originalpckg./Stück. Stand. pack/quantity
24 205 01	80 x 80	250	10
24 205 03	100 x 100	500	10
24 205 05	120 x 120	1 000	1
24 205 10	150 x 150	2 000	1
24 205 21	210 x 210	6 000	1
24 205 32	260 x 260	12 000	1


**Zylinder**  
mit Knopfdeckel, Rand  
poliert

**Bocaux**  
avec couvercle à  
bouton, bord poli

**Cylinders**  
with knobbed lid,  
polished rim

**Tarros cilíndricos**  
tapa con bola, bordes  
pulidos

**NEU/NEW/NOUVEAU/NUEVO**



Best.-Nr. Cat. No.	Inhalt Capacity ml	ID ID mm	AD OD mm	h mm	Originalpckg./Stück. Stand. pack/quantity
10 088 06	1 000	100	107	195	1
10 088 07	2 000	100	107	315	1

**DURAN® safe-**  
Präparatengläser

**DURAN® safe-**  
Verres à préparation

**DURAN® safe-**  
Specimen jars

**DURAN® safe-**  
Tarros cilíndricos para  
preparaciones

Best.-Nr. Cat. No.	Inhalt Capacity ml	ID ID mm	AD OD mm	h mm	Originalpckg./Stück. Stand. pack/quantity
DURAN® safe-Glasgefäß/DURAN® safe glass vessel					
24 203 5002	1 000	100	107	180	1
10 088 05	2 000	100	107	300	1
DURAN® safe-Verschluss/DURAN® safe closure					
10 088 04		100		46	1