


Best.-Nr. Cat. No.	h + h <sub>1</sub> mm	h x d mm	NS	Originalpckg./Stück Stand. pack/quantity
24 465 59	300	250 x 185	34/35	1
24 465 61	350	300 x 215	34/35	1
24 465 69	550	500 x 315	34/35	1




**Planflansch-Glocke**  
mit Bohrung im Hals,  
vakuumfest

**Bell jars**  
with neck bore, for  
vacuum use

**Cloche à bride plane**  
avec tubulure, pour  
emploi sous vide

**Campana con brida  
plana**  
con tubuladura,  
resistentes al vacío

Best.-Nr. Cat. No.	Inhalt Capacity ml	größter d max. d mm	h mm	Originalpckg./Stück Stand. pack/quantity
21 401 54	1000	120	465	10



**DURAN®-Sedimentier-  
gefäß nach Imhoff**  
graduirt, DIN 12 672,  
Teil 1

**Sedimentation cone,  
made of DURAN®,  
Imhoff type**  
graduated, DIN 12 672,  
Part 1

**Cône à sédimentation  
en DURAN® d'après  
Imhoff**  
gradué, DIN 12 672,  
partie 1

**Conos de sedimentación  
según Imhoff, en vidrio  
DURAN®**  
graduados, DIN 12 672,  
Parte 1