


KALK-SODA/SODA-LIME/  
SODO-CALCIQUE/CAL-SODA

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN thread GL	d mm	h mm	Originalpckg./Stück Stand. pack/quantity
Klarglas/clear glass/verre clair/vidrio transparente					
23 810 24 5	100	32 <sup>1</sup>	49	119	30
23 810 36 5	250	32 <sup>1</sup>	64	155	10
23 810 44 5	500	32 <sup>1</sup>	77	186	10
23 810 54 5	1 000	45 <sup>1</sup>	96	223	10
Braunglas/amber glass/verre brun/vidrio topacio					
23 816 24 5	100	32 <sup>1</sup>	49	119	10
23 816 36 5	250	32 <sup>1</sup>	64	155	10
23 816 44 5	500	32 <sup>1</sup>	77	186	10
23 816 54 5	1 000	45 <sup>1</sup>	96	223	10

<sup>1</sup> Hohe Form/highform/forme haute/forma alta

**Vierkant-Schraubflaschen**, enghalsig mit DIN-Gewinde, aus Kalk-Soda-Glas, **komplett** mit Ausgießring und Staubschutzkappe

**Square bottles**, narrow neck, with DIN thread, made of **soda-lime glass**, **complete** with dust cap and pouring ring



**Flacons carrés**, à col étroit, à filetage DIN, en **verre sodo-calcique**, **complets** avec bague de déversement et capuchon de protection anti-poussière


**Frasco cuadrangular**, cuello estrecho, con rosca DIN, en **vidrio cal-soda**, **completo** con anillo de vertido y tapa de protección contra el polvo

KALK-SODA/SODA-LIME/  
SODO-CALCIQUE/CAL-SODA

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN thread GL	d mm	h mm	Originalpckg./Stück Stand. pack/quantity
Klarglas/clear glass/verre clair/vidrio transparente					
23 820 24	100	32 <sup>1</sup>	49	109	30
23 820 36	250	45 <sup>1</sup>	64	146	10
23 820 44	500	54 <sup>1, 2</sup>	77	173	10
23 820 54	1 000	60 <sup>1, 2</sup>	97	213	10
Braunglas/amber glass/verre brun/vidrio topacio					
23 826 24	100	32 <sup>1</sup>	49	109	10
23 826 36	250	45 <sup>1</sup>	64	146	10
23 826 44	500	54 <sup>1, 2</sup>	77	173	10
23 826 54	1 000	54 <sup>1, 2</sup>	97	213	10


<sup>1</sup> Kurze Form/short form/forme courte/forma baja

Passende Kappe siehe Seite 34, 35.  
For spare caps see page 34, 35.  
Capuchons, voir page 34, 35.  
Tapas adecuadas: ver pág. 34, 35.

<sup>2</sup>Nicht nach DIN<sup>2</sup>Non DIN-size<sup>2</sup>Non conforme aux dimensions DIN<sup>2</sup>No según DIN

**Vierkant-Schraubflaschen**, weithalsig mit DIN-Gewinde, aus Kalk-Soda-Glas, **ohne** Ausgießring und Staubschutzkappe

**Square bottles**, wide neck, with DIN thread, made of **soda-lime glass**, **without** cap and pouring ring



**Flacons carrés**, à col large, à filetage DIN, en **verre sodo-calcique**, **sans** bague de déversement ni capuchon de protection anti-poussière

**Frasco cuadrangular**, cuello ancho, con rosca DIN, en **vidrio cal-soda**, **sin** anillo de vertido y tapa de protección contra el polvo