


Best.-Nr. Cat. No.	DN	h mm	ID Flansch ID flange mm	AD Flansch OD flange mm	ca. Vol. approx. vol. l	Tubus Tubulature	Originalpckg./Stück Stand. pack/quantity
Klarglas/clear glass							
24 782 46	100	174	118,5	153 +/-2	0,7	24/29	1
24 782 57	150	239	172	215 +/-2	2,4	24/29	1
24 782 61	200	296	222	270 +/-2	5,8	24/29	1
24 782 66	250	344	272	320 +/-2	10,5	24/29	1
24 782 69	300	420	332	380 +/-2	18,5	24/29	1



**Vakuum-Exsikkatoren**  
mit NOVUS NS-Tubus  
im Deckel, mit Hahn,  
mit Planflansch,  
DIN 12 491


**Desiccators**  
with flat flange and  
NOVUS NS junction.  
Outlet in lid, with  
stopcock, DIN 12 491



**Desiccateurs à vide**  
avec tubulure RIN  
NOVUS au couvercle,  
à rodage plan, avec  
robinet, DIN 12 491


**Desecadores al vacío**  
con tubuladura NS  
NOVUS en la tapa, con  
brida plana, con llave,  
DIN 12 491

Best.-Nr. Cat. No.	DN	h mm	ID Flansch ID flange mm	AD Flansch OD flange mm	ca. Vol. approx. vol. l	Originalpckg./Stück Stand. pack/quantity
Klarglas/clear glass						
24 781 46	100	187	118,5	153 +/-2	0,7	1
24 781 57	150	252	172	215 +/-2	2,4	1
24 781 61	200	309	222	270 +/-2	5,8	1
24 781 66	250	357	272	320 +/-2	10,5	1
24 781 69	300	433	332	380 +/-2	18,5	1



**Vakuum-Exsikkatoren**  
mit Planflansch, ohne  
Anschluss, mit Knopf-  
deckel

**Desiccators**  
with plane flange, with-  
out junction, with knob-  
bed lid



**Desiccateurs à vide**  
à rodage plan, sans  
raccord, avec couvercle  
à bouton

**Desecadores al vacío**  
con brida plana, sin co-  
nexión, con tapa a bola