


Best.-Nr. Cat. No.	DN	h	h <sub>1</sub>	ID Flansch ID flange	AD Flansch OD flange	ca. Vol. approx. vol.	Original- pckg./Stück Stand. pack/quantity
		mm	mm	mm	mm	l	
Klarglas/clear glass							
24 772 57	150	154	81	172	215 +/-2	2,4	1
24 772 61	200	202	115	222	270 +/-2	5,8	1
24 772 66	250	235	120	272	320 +/-2	10,5	1
24 772 69	300	283	150	332	380 +/-2	18,5	1

**Exsikkator-Unterteile**  
mit Planflansch,  
Gewindeanschluss,  
Typ MOBILEX

**Parties inférieures  
de dessiccateurs**  
à rodage plan, raccord à  
filetage, type MOBILEX

**Desiccator bases**  
with flat flange,  
screw thread outlet,  
MOBILEX type

**Base para desecadores**  
con brida plana, con  
tubuladura a rosca, tipo  
MOBILEX




Best.-Nr. Cat. No.	DN	h	ID Flansch ID flange	AD Flansch OD flange	Originalpckg./ Stück Stand. pack/ quantity
		mm	mm	mm	
Klarglas/clear glass					
24 410 46	100	75	118,5	153 +/-2	1
24 410 57	150	98	172	215 +/-2	1
24 410 61	200	107	222	270 +/-2	1
24 410 66	250	122	272	320 +/-2	1
24 410 69	300	150	332	380 +/-2	1

**Exsikkator-Deckel**  
mit Knopf, passend  
für alle Unterteile

**Couvercles pour  
dessiccateurs**  
à bouton, adaptables  
à toutes les parties  
inférieures

**Desiccator lids**  
with knob, fits all  
bases

**Tapa para desecador**  
con bola, adecuado  
para todas las bases



Best.-Nr. Cat. No.	DN	h	ID Flansch ID flange	AD Flansch OD flange	NS	Originalpckg./ Stück Stand. pack/ quantity
		mm	mm	mm		
Klarglas/clear glass						
24 420 46	100	62	118,5	153 +/-2	24/29	1
24 420 57	150	85	172	215 +/-2	24/29	1
24 420 61	200	94	222	270 +/-2	24/29	1
24 420 66	250	109	272	320 +/-2	24/29	1
24 420 69	300	137	332	380 +/-2	24/29	1

**Exsikkator-Deckel**  
für normale NS-Hähne,  
Typ NOVUS,  
DIN 12 491

**Couvercles pour  
dessiccateurs**  
pour robinets RIN  
normaux type NOVUS,  
DIN 12 491

**Desiccator lids**  
outlet for normal NS  
stopcocks NOVUS type,  
DIN 12 491

**Tapa para desecador**  
para llaves normales  
NS tipo NOVUS,  
DIN 12 491