

**DURAN®-Zentrifugen-
gläser mit Rundboden**
DIN 58 970, Teil 2


**Tubes à centrifuger
à fond rond en verre
DURAN®**
DIN 58 970, section 2

**DURAN® centrifuge
tubes, round bottom**
DIN 58 970, part 2

**DURAN®-Tubos para
centrífuga con fondo
redondo**
DIN 58 970, parte 2

Best.-Nr. Cat. No.	d x h mm	Nenninhalt Nominal capacity ml	Originalpckg./Stück Stand. pack/quantity
21 601 10	12 x 100	6	50
21 601 11	16 x 100	12	50
21 601 14	24 x 100	25	10
21 601 17	34 x 100	50	10
21 601 24 ¹	40 x 115	80	10
21 601 26	44 x 100	80	10
21 601 36 ¹	56 x 147	250	10

¹ Nicht nach DIN/Non-DIN size/Non conforme aux dimensions DIN/No según DIN




**DURAN®-Zentrifugen-
glas mit Spitzboden**
Winkel 30°

**Tube à centrifuger
cylindroconique
en verre DURAN®**
angle 30°

**DURAN® centrifuge
tube, conical**
angle 30°

**DURAN®-Tubo para
centrífuga con fondo
cónico, ángulo 30°**

Best.-Nr. Cat. No.	d x h mm	Nenninhalt Nominal capacity ml	Originalpckg./Stück Stand. pack/quantity
24 263 09	16 x 100	12	50



**DURAN®-Zentrifugen-
glas mit Spitzboden**
Winkel 60°

**Tubes à centrifuger
cylindroconiques en
verre DURAN®**
angle 60°

**DURAN® centrifuge
tube, conical**
angle 60°

**DURAN®-Tubos para
centrífuga con fondo
cónico, ángulo 60°**

Best.-Nr. Cat. No.	d x h mm	Nenninhalt Nominal capacity ml	Originalpckg./Stück Stand. pack/quantity
21 611 14	24 x 100	25	10
21 611 17	34 x 100	50	10