


Best.-Nr. Cat. No.	Inhalt Capacity ml	l mm	b mm	h mm	Originalpckg./Stück Stand. pack/quantity
21 571 58	1 200	275	132	56	10

Kulturkolben nach Roux
Hals konisch, exzentrisch

Ballon à culture selon Roux
col conique, excentrique

Culture flask, Roux-type
conical neck, excentric

Matraz para cultivos según Roux
cuello cónico excéntrico




Best.-Nr. Cat. No.	Inhalt Capacity ml	l mm	b mm	h mm	d mm	Originalpckg./Stück Stand. pack/quantity
21 551 71	4 000	370	200	90	50	1

Penicillinkolben

Boîte à pénicilline à bord droit

Penicillin flask

Matraz para penicilina



Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d ₁ mm	h mm	Best.-Nr. für pass. Kapsenberg-Kappe Cat.-Nr. for matching Kapsenberg caps	Originalpckg./Stück Stand. pack/quantity
21 421 17	50	40	18	107	29 010 11	10
21 421 24	100	40	18	150	29 010 11	10
21 421 32 ¹	200	50	18	175	29 010 11	10

¹ DIN 38 411, Teil 6
¹ DIN 38 411, part 6
¹ DIN 38 411, partie 6
¹ DIN 38 411, apartado 6

Kulturflaschen
Rand gerade
für Kapsenberg-Kappen

Flacons à culture à bord droit
pour emploi avec capuchons
de Kapsenberg

Culture bottles
with straight rim
for Kapsenberg caps

Frascos para cultivos,
borde recto
para casquetes de
Kapsenberg