


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Originalpckg./Stück Stand. pack/quantity
21 771 24	100	64	38	114	10
21 771 32	200	79	38	138	10
21 771 36	250	85	38	149	10
21 771 39	300	87	38	161	10
21 771 44	500	105	38	183	10
21 771 54	1 000	131	38	232	10
21 771 63	2 000	166	38	305	10




**Kulturkolben, Erlenmeyerform, Hals gerade für Metall-Kappen**

**Flacons pour milieux de culture, col droit pour capuchons en métal**

**Culture flasks, Erlenmeyer shape, straight neck for metal caps**

**Frascos para cultivos, forma Erlenmeyer, cuello recto para casquetes de metal**

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Originalpckg./Stück Stand. pack/quantity
21 773 24	100	50	38	125	10
21 773 39	300	71	38	170	10
21 773 44	500	83	38	208	10
21 773 54	1 000	105	38	240	10




**Nährbodenflaschen Hals gerade für Metall-Kappen**

**Flacons pour milieux de culture, col droit pour capuchons en métal**

**Culture media bottles straight neck for metal caps**

**Frascos para cultivos, cuello recto para casquetes de metal**

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Originalpckg./Stück Stand. pack/quantity
21 774 62	1 800	200	38	175	2



**Kulturkolben nach Fernbach, konische Form, Hals gerade für Metall-Kappen**

**Ballon à culture selon Fernbach, forme conique, col droit pour capuchons en métal**

**Culture flask, Fernbach type, conical shape straight neck for metal caps**

**Matraces para cultivos, según Fernbach, forma cónica, cuello recto para capuchones de metal**