


AR-GLAS®/AR-GLASS®/
VERRE AR®/VIDRIO AR®

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Farbkenn- zeichnung Colour code DIN 12 621	Ablaufzeit Run-out time s	Gesamt- länge Overall length mm	Original- pckg./Stück Stand. pack/ quantity
24 338 01 ¹	1	0,007	blau blue	5– 9	270	12
24 338 02	2	0,01	orange	5– 9	330	12
24 338 07	5	0,015	weiß white	7–11	400	12
24 338 08	10	0,02	rot red	8–12	440	12
24 338 12	20	0,03	gelb yellow	9–13	510	6
24 338 14	25	0,03	blau blue	10–15	520	6
24 338 17	50	0,05	rot red	13–18	540	6
24 338 24	100	0,08	gelb yellow	25–30	585	6

¹ Ohne Flüssigkeitsreservoir¹ Without liquid reservoir¹ Sans réservoir de liquide¹ Sin depósito reserva de líquido**AR-GLAS®-Vollpipetten**Klasse AS, DIN 12 691,
ISO 648, Wartezeit 15 s**Pipettes jaugées**classe AS, DIN 12 691,
ISO 648, temps
d'attente 15 secondes**AR-GLASS® bulb****pipettes, graduated**
class AS, DIN 12 691,
ISO 648, waiting time
15 seconds**Pipetas aforadas**clase AS, DIN 12 691,
ISO 648, tiempo de
espera 15 segundos