

Best.-Nr. Cat. No. DURAN®	Best.-Nr. Cat. No. Kalk-Soda Soda-lime	d mm	h mm	~ h <sub>1</sub> mm	d <sub>1</sub> mm	Pass. Papierfilter Match. filter paper Ø mm	Originalpckg./Stück Stand. pack/quantity
21 351 23	23 351 23	35	60	35	6	45– 55	10
21 351 28	23 351 28 <sup>2</sup>	45	80	45	6	55– 70	10
21 351 33	23 351 33	55	95	55	8	70– 90	10
21 351 38 <sup>2</sup>	23 351 38 <sup>2</sup>	70	120	70	8	110–125	10
21 351 41 <sup>2</sup>	23 351 41 <sup>2</sup>	80	140	80	10	125–150	10
21 351 46	23 351 46	100	180	100	10	150–185	10
21 351 51 <sup>1,2</sup>	–	120	210	120	16	185–240	10
21 351 57	–	150	250	150	16	240–270	10
21 351 59 <sup>1,2</sup>	–	180	290	150	20	270–320	1
21 351 61 <sup>2</sup>	–	200	325	175	26	320–385	1
21 351 66 <sup>1,2</sup>	–	250	360	175	30	385–400	1
21 351 69 <sup>2</sup>	–	300	409	175	30	500	1

<sup>1</sup> Nicht nach DIN

<sup>1</sup> Non-DIN size

<sup>1</sup> Non conforme aux dimensions DIN


<sup>1</sup> No según DIN

<sup>2</sup> Nicht nach ISO

<sup>2</sup> Non-ISO size

<sup>2</sup> Non conforme aux dimensions ISO

<sup>2</sup> No según ISO



**Trichter**  
mit kurzem Stiel,  
DIN 12 445, ISO 4798




**Entonnoir**  
à tige courte,  
DIN 12 445, ISO 4798

**Funnels**  
with short stem,  
DIN 12 445, ISO 4798

**Embudos**  
rama corta,  
DIN 12 445, ISO 4798

Best.-Nr. Cat. No. DURAN®	Best.-Nr. Cat. No. Kalk-Soda Soda-lime	d mm	h mm	~ h <sub>1</sub> mm	d <sub>1</sub> mm	Originalpckg./Stück Stand. pack/quantity
21 354 33	23 354 33	55	60	30	20	10
21 354 38	23 354 38	70	70	30	22	10
21 354 41	23 354 41	80	80	30	24	10
21 354 46	23 354 46	100	100	30	26	10
21 354 51	–	120	105	30	35	10
21 354 55	–	160	140	30	35	1
21 354 61	–	200	170	30	40	1



**Pulvertrichter**  
mit kurzem, weitem Stiel



**Entonnoir à poudre**  
à tige courte et large

**Powder funnels**  
with short, wide stem

**Embudos para polvo**  
rama corta y ancha

Best.-Nr. Cat. No. DURAN®	Best.-Nr. Cat. No. Kalk-Soda Soda-lime	d mm	h mm	~ h <sub>1</sub> mm	d <sub>1</sub> mm	Pass. Papierfilter Match. filter paper Ø mm	Originalpckg./Stück Stand. pack/quantity
21 353 33	23 353 33 <sup>1</sup>	55	190	150	6 <sup>3</sup>	70– 90	10
21 353 38 <sup>2</sup>	23 353 38 <sup>1,2</sup>	70	200	150	6 <sup>3</sup>	110–125	10
21 353 41 <sup>2</sup>	23 353 41 <sup>1,2</sup>	80	210	150	6 <sup>3</sup>	125–150	10
–	23 353 46 <sup>1</sup>	100	240	150	9	150–185	10

<sup>1</sup> Nicht nach DIN

<sup>1</sup> Non-DIN size

<sup>1</sup> Non conforme aux dimensions DIN

<sup>1</sup> No según DIN

<sup>2</sup> Nicht nach ISO

<sup>2</sup> Non-ISO size

<sup>2</sup> Non conforme aux dimensions ISO


<sup>2</sup> No según ISO

<sup>3</sup>Anhaben gelten nur für DURAN®-Ausführung

<sup>3</sup>Measurements only valid for DURAN®-type

<sup>3</sup>Les références indiquées ci-dessus ne se réfèrent qu'aux entonnoirs faits en DURAN®

<sup>3</sup>Datos son validos únicamente para DURAN®



**Trichter (Bunsenrichter)**  
mit langem Stiel,  
DIN 12 446, ISO 4798



**Entonnoir**  
à tige longue,  
DIN 12 446, ISO 4798

**Funnels**  
with long stem,  
DIN 12 446, ISO 4798

**Embudos**  
rama longa,  
DIN 12 446, ISO 4798