


Analysentrichter
für Schnellfiltration

Entonnoir d'analyse
pour filtration rapide

Analytical funnels
for rapid filtration

Embudos de análisis
para filtración rápida

Best.-Nr. Cat. No.	d mm	h mm	h ₁ mm	d ₁ mm	Pass. Papier-Rundfilter Match. filter paper Ø mm	Originalpckg./Stück Stand. pack/quantity
21 331 37	65	200	150	9	70– 90	10
21 331 41	80	210	150	9	110–125	10
21 331 48	110	265	180	9	150–185	10




Trichter gerippt

Entonnoir à rainures

Funnels ribbed

Embudos estriados

Best.-Nr. Cat. No.	d mm	h mm	h ₁ mm	d ₁ mm	Pass. Papier-Rundfilter Match. filter paper Ø mm	Originalpckg./Stück Stand. pack/quantity
21 352 38	70	125	70	8	110–125	10
21 352 41	80	140	80	10	125–150	10
21 352 46	100	180	100	10	150–185	10
21 352 57	150	266	150	16	240–270	10
21 352 61	200	326	175	26	320–385	1




Filtertrichter
Kegelform

Entonnoirs filtrants
coniques

Filter funnels
conical shape

Embudos filtrantes
forma cónica

Best.-Nr. Cat. No.	Kennzeichnung Identification mark	Porosität Porosity	Inhalt Capacity ml	d mm	d ₁ mm	d ₂ mm	h mm	Originalpckg./Stück Stand. pack/quantity
25 854 03	39D	3	25	55	25	8	100	1
25 854 04	39D	4	25	55	25	8	100	1



Filternutschen

Entonnoirs filtrants

Filter funnels

Embudos filtrantes

Best.-Nr. Cat. No.	Kennzeichnung Identification mark	Porosität Porosity	Inhalt Capacity ml	Platten-Ø mm	d ₁ mm	h mm	ID ID	Originalpckg./Stück Stand. pack/quantity
25 852 0X ¹	3D	1–5	50	35	10	130	40	1
25 852 1X ¹	11D	1–5	75	45	10	130	56	1
25 852 2X ¹	17D	1–5	125	60	10	140	71	1
25 852 3X ¹	25D	1–5	500	95	22	200	106	1
25 852 4X ¹	26D	1–5	1 000	120	22	320	135	1
25 852 6X ¹	151D	1–4	4 000	175	30	400	202	1

¹ In die Best.-Nr. ist statt X die gewünschte Porosität einzutragen, siehe Seite 147, 148.

¹ The required porosity should be inserted in place of X in the catalogue number, see page 169, 170.

¹ Dans le numéro d'article, remplacer X par la porosité voulue, voir page 191, 192.

¹ En el número de referencia se debe poner en lugar de la X la porosidad deseada, ver página 213, 214.