


Best.-Nr. Cat. No.	Inhalt Capacity ml	Hälsa Neck NS	d mm	Originalpckg./Stück Stand. pack/quantity
24 709 44	500	19/26	87	1
24 709 54	1 000	24/29	113	1
24 709 63	2 000	29/32	135	1
24 709 73	5 000	34/35	185	1




**DURAN®-Woulff'sche
Flaschen**
mit 3 Hälsen NS,
DIN 12 480, vakuumfest

**Flacons de Woulff
en DURAN®**
à 3 cols RIN,
DIN 12 480,
pour emploi sous vide

DURAN®-Woulff bottles
with 3 necks,
DIN 12 480,
for vacuum use

**Frascos de Woulff
en vidrio DURAN®**
con 3 bocas NS,
DIN 12 480, resistentes
al vacío

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hälsa Neck NS	Boden Bottom NS	d mm	Originalpckg./Stück Stand. pack/quantity
24 710 44	500	19/26	19/26	87	1
24 710 54	1 000	24/29	19/26	113	1
24 710 63	2 000	29/32	19/26	135	1
24 710 73	5 000	34/35	29/32	185	1



**DURAN®-Woulff'sche
Flaschen**
mit 3 Hälsen NS und
Bodentubus,
DIN 12 481, vakuumfest

**Flacons de Woulff
en DURAN®**
à 3 cols RIN et tubulure
de fond, DIN 12 481,
pour emploi sous vide

DURAN®-Woulff bottles
with 3 ground necks NS
and bottom tubulature,
DIN 12 481, for vacuum
use

**Frascos de Woulff
en vidrio DURAN®**
con 3 bocas NS y
tubuladura, DIN 12 481,
resistentes al vacío