


| Best.-Nr. Cat. No. | Inhalt Capacity ml | NS Socket size | d mm | h mm | Originalpkg./Stück Stand. pack/quantity |
|------------------------|--------------------------|-------------------|---------|---------|--|
| 24 193 13 | 25 | 14/23 | 42 | 75 | 10 |
| 24 193 20 | 50 | 14/23 | 51 | 85 | 10 |
| 24 193 17 | 50 | 19/26 | 51 | 85 | 10 |
| 24 193 18 ¹ | 50 | 24/29 | 51 | 85 | 10 |
| 24 193 19 | 50 | 29/32 | 51 | 85 | 10 |
| 24 193 24 | 100 | 19/26 | 64 | 105 | 10 |
| 24 193 26 ¹ | 100 | 24/29 | 64 | 105 | 10 |
| 24 193 27 | 100 | 29/32 | 64 | 105 | 10 |
| 24 193 32 ¹ | 200 | 29/32 | 79 | 131 | 10 |
| 24 193 36 ¹ | 250 | 24/29 | 85 | 140 | 10 |
| 24 193 37 | 250 | 29/32 | 85 | 140 | 10 |
| 24 193 38 | 250 | 45/40 | 85 | 140 | 10 |
| 24 193 39 ¹ | 300 | 29/32 | 87 | 156 | 10 |
| 24 193 44 ¹ | 500 | 24/29 | 105 | 175 | 10 |
| 24 193 46 | 500 | 29/32 | 105 | 175 | 10 |
| 24 193 47 | 500 | 45/40 | 105 | 175 | 10 |
| 24 193 54 ¹ | 1 000 | 24/29 | 131 | 220 | 10 |
| 24 193 56 | 1 000 | 29/32 | 131 | 220 | 10 |
| 24 193 57 | 1 000 | 45/40 | 131 | 220 | 10 |

¹ Nicht nach DIN¹ Non-DIN size¹ Taille spéciale, non conforme aux dimensions DIN¹ Tamaño especial, no según DIN

Erlenmeyer-Kolben
mit Normschliff,
DIN 12 387


Flasks, Erlenmeyer
with standard ground
socket, DIN 12 387



Fioles Erlenmeyer
à joint rodé standard,
DIN 12 387

Matraces Erlenmeyer
con esmerilado
normalizado, DIN 12 387

| Best.-Nr. Cat. No. | Inhalt Capacity ml | NS Socket size | d mm | h mm | Originalpkg./Stück Stand. pack/quantity |
|------------------------|--------------------------|-------------------|---------|---------|--|
| 24 171 19 | 50 | 29/32 | 51 | 90 | 10 |
| 24 171 24 | 100 | 19/26 | 64 | 110 | 10 |
| 24 171 26 ¹ | 100 | 24/29 | 64 | 110 | 10 |
| 24 171 27 | 100 | 29/32 | 64 | 110 | 10 |
| 24 171 36 ¹ | 250 | 24/29 | 85 | 140 | 10 |
| 24 171 37 | 250 | 29/32 | 85 | 140 | 10 |
| 24 171 44 ¹ | 500 | 24/29 | 105 | 170 | 10 |
| 24 171 46 | 500 | 29/32 | 105 | 170 | 10 |
| 24 171 54 ¹ | 1 000 | 24/29 | 131 | 200 | 10 |
| 24 171 56 | 1 000 | 29/32 | 131 | 200 | 10 |
| 24 171 63 | 2 000 | 29/32 | 166 | 250 | 10 |

¹ Nicht nach DIN¹ Non-DIN size¹ Taille spéciale, non conforme aux dimensions DIN¹ Tamaño especial, no según DIN

Stehkolben
mit Normschliff,
DIN 12 348


Flasks, flat bottom
with standard ground
socket, DIN 12 348



Ballons à fond plat,
longueur moyenne
avec joint rodé standard,
DIN 12 348

Matraces fondo plano,
largo mediano
con esmerilado normali-
zado, DIN 12 348

| Best.-Nr. Cat. No. | Inhalt Capacity ml | NS Socket size | d mm | h mm | Originalpkg./Stück Stand. pack/quantity |
|-----------------------|--------------------------|-------------------|---------|---------|--|
| 24 120 27 | 100 | 29/32 | 60 | 110 | 10 |
| 24 120 37 | 250 | 29/32 | 81 | 140 | 10 |
| 24 120 46 | 500 | 29/32 | 101 | 170 | 10 |
| 24 120 56 | 1 000 | 29/32 | 126 | 210 | 10 |



Verdampferkolben
mit Normschliff,
birnenförmig

Flasks, evaporating
with standard ground
socket, pear shape



Ballons d'évaporation
avec joint rodé
standard, piriformes

Matraces para evaporar
con esmerilado
normalizado, forma de
pera