





Best.-Nr. Cat. No.	Kennzeichnung Designation	d mm	h mm	Kern NS Cone size	Originalpckg./Stück Stand. pack/quantity
24 523 55	HP 10	10	75	24/29	1
24 528 56	HQ 10	10	75	29/32	1
24 535 49	HS 10	10	65	24/29	1
24 540 51	HT 10	10	65	29/32	1
24 541 54	HT 16	16	85	29/32	1





KPG-Lagerhülsen
austauschbar, Lauffläche
geschliffen und poliert,
mit Normschliffkern

**Joints d'agitateur à
palier cylindrique**
interchangeables, surface
de palier rodée et polie,
avec joint mâle à rodage
standard

KPG stirrer bearings
interchangeable, bearing
surface ground and
polished, with standard
ground cone

Cierres de agitación KPG
intercambiables, superficie
de rozamiento pulida,
con esmerilado macho
normalizado

Best.-Nr. Cat. No.	Kennzeichnung Designation	d mm	h mm	h ₁ mm	Originalpckg./Stück Stand. pack/quantity
24 500 42	HA 10	10	80	65	1
24 505 51	HB 10	10	120	75	1
24 506 57	HB 16	16	150	90	1



KPG-Lagerhülsen
austauschbar, Lauffläche
geschliffen und poliert

**Joints d'agitateur
à palier cylindrique**
interchangeables, surface
de palier rodée et polie

KPG stirrer bearings
interchangeable,
bearing surface ground
and polished

Cierres de agitación KPG
intercambiables, superficie
de rozamiento pulida