



Best.-Nr. Cat. No.	Kenn- zeichnung Designation	d mm	h mm	Din-Gewinde DIN thread GL	Kern NS Cone size	Originalpckg./Stück Stand. pack/quantity
24 750 08	HB 10	10	75	32	24/29	1
24 750 09	HB 10	10	75	32	29/32	1

KPG-Lagerhülsen
austauschbar, Lauffläche
geschliffen und poliert,
mit Gewinderohr und
Normschliffkern

**Joints d'agitateur à
palier cylindrique**
interchangeables, surface
de palier rodée et polie,
avec tube fileté et joint
mâle à rodage standard

KPG stirrer bearings
interchangeable, bearing
surface ground and
polished, with screw-
thread and standard
ground cone

Cierres de agitación KPG
intercambiables,
superficie de rozamiento
pulida, provistos de rosca
y esmerilado macho
normalizado



Best.-Nr. Cat. No.	Kenn- zeichnung Designation	Welle		Gesamt- länge Overall length mm	Passend für Kolbenhals Suitable for flask neck ID mm	Original- pckg./Stück Stand. pack/ quantity
		d mm	l mm			
24 573 74	WG 10	10	160	320	60	1
24 573 77	WG 10	10	160	370	60	1
24 573 84	WG 10	10	160	410	60	1
24 573 86	WG 10	10	160	440	60	1
24 578 74	WH 10	10	160	320	60	1
24 578 77	WH 10	10	160	370	60	1
24 578 84	WH 10	10	160	410	60	1
24 578 86	WH 10	10	160	440	60	1
24 583 74	WS 10	10	160	320	25	1
24 583 77	WS 10	10	160	370	25	1
24 583 84	WS 10	10	160	410	25	1
24 583 86	WS 10	10	160	440	25	1

**KPG-Rührerwellen
Ø 10 mm**
austauschbar, Lauffläche
geschliffen und poliert

**Tiges calibrées pour
agitateurs en verre,
Ø 10 mm**
interchangeables, surface
de palier rodée et polie

**KPG stirrer shafts,
Ø 10 mm**
interchangeable,
bearing surface ground
and polished

**Agitadores KPG
Ø 10 mm**
intercambiables, superficie
de rozamiento esmerilada
y pulida