



DURAN® SAFETY JOINTS

Securing joints by using our KECK™-products is a well-known and proven method.

In order to meet high safety demands, e.g. when dealing with hazardous substances or pressure differences, DURAN® Safety Joints offer very good protection against accidental loosening of the connection. Furthermore, DURAN® Safety Joints offer a system to separate jammed cones and sockets. This minimizes the risk of breakage and related threats.

APPLICATIONS

- Handling hazardous substances
- Avoiding the risk of stuck joints when working grease-free
- For mechanically stressed connections, e.g. as a result of the assembly situation or vibration
- For tight and secure closure

FEATURES

- Borosilicate glass 3.3
- Reliable connection
- Sealing properties as good as DURAN® standard joints
- Easy and safe separation
- Compatible with DURAN® cones and sockets according to DIN 12 249

ORDER INFORMATION

DURAN® THREADED SOCKET

Nominal size	Pack unit	Cat.-No.
NS 14/23	10	28 659 14 04
NS 19/26	10	28 659 19 01
NS 24/29	10	28 659 24 09
NS 29/32	10	28 659 29 06
NS 34/35	10	28 659 34 05
NS 45/40	10	28 659 45 04



DURAN® CONE WITH RIM

Nominal size	Pack unit	Cat.-No.
NS 14/23	10	28 618 14 02
NS 19/26	10	28 618 19 08
NS 24/29	10	28 618 24 07
NS 29/32	10	28 618 29 04
NS 34/35	10	28 618 34 03
NS 45/40	10	28 618 45 02



DURAN® CONE WITH TIP AND RIM

Nominal size	Pack unit	Cat.-No.
NS 19/26	10	28 619 19 09
NS 29/32	10	28 619 29 05



ACCESSORIES FOR DURAN® SAFETY JOINTS AND FOR RODAVISS®*

SCREWCAPS WITH APERTURE

Open screwcaps ensure the connection of the cone and the socket.
By loosening the connection, the parts are automatically separated.

For nominal size	Pack unit	Cat. No.
NS 14	10	28 621 50 01
NS 19	10	28 621 52 07
NS 24	10	28 621 53 01
NS 29	10	28 621 54 04
NS 34	10	28 621 55 07
NS 45	10	28 621 57 04



O-RING

O-Rings (Nitrile) for DURAN® Safety Joints enable the connection
of cone and socket along with the screw cap.

For nominal size	Pack unit	Cat. No.
NS 14	10	28 621 68 03
NS 19	10	28 621 70 02
NS 24	10	28 621 72 08
NS 29	10	28 621 74 05
NS 34	10	28 621 75 08
NS 45	10	28 621 77 05



LOOSENING RING

Loosening Rings (polyamide) have very good sliding properties.
By distributing force evenly, they allow stuck joints to be safely disconnected.


For nominal size	Pack unit	Cat. No.
NS 14	10	28 621 64 09
NS 19	10	28 621 69 06
NS 24	10	28 621 71 05
NS 29	10	28 621 66 06
NS 34	10	28 621 73 02
NS 45	10	28 621 67 09



SCREWCAPS CLOSED

Closed screw caps enable apparatus to be closed securely after use.
PTFE-coated silicon seals provide maximum chemical resistance.

For nominal size	Pack unit	Cat. No.
NS 14	10	28 621 42 02
NS 19	10	28 621 43 05
NS 24	10	28 621 44 08
NS 29	10	28 621 45 02
NS 34	10	28 621 46 05
NS 45	10	28 621 47 08



* RODAVISS® is a registered trademark of SAS SOUFFLAGE ARTISANAL DU VERRE



DWK Life Sciences GmbH
Otto-Schott-Str. 21
97877 Wertheim / Main
Germany

Phone: +49 9342 - 802 214
Fax: +49 9342 - 802 250
glassblowers@DWK-LifeSciences.com
www.DWK-LifeSciences.com

