



**DURAN  
WHEATON  
KIMBLE**

Excellence in your hands

# WHEATON®

## AUTOMATIC CRIMPING TOOL BATTERY POWERED



### Battery Powered Crimping Tools

Our battery powered crimpers and decappers fit standard laboratory vials and seals with aluminum or two-part aluminum/ steel caps. The electronically adjustable cycle is accurate and reproducible and can be adjusted with the control pad on top of the unit. The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools.

Rechargeable: The tools can be recharged in one to two hours with the universal power supply included. The power supply includes a plug set for operation in most countries.

Product Highlights	
• Powerful motor with a fast cycle	• Intuitive LCD screen makes adjusting the crimp easy
• Long lasting battery with Lithium Ion technology - up to 800 crimps per charge	• Advanced charging system including warning icon when recharging is necessary
• Can be operated while charging	• Optional base and mounting kit
• For flip-off caps, available in 13 mm and 20 mm sizes*	

\*Not for use with Flip-up caps

# AUTOMATIC CRIMPING TOOL — BATTERY POWERED

## Battery Powered Crimping Tool

Cat. No.	Description	Size (mm)	Qty
W225808	Battery Powered Crimper	8	1
W225811	Battery Powered Crimper	11	1
W225812	Battery Powered Decapper	11	1
W225813	Battery Powered Crimper	13	1
W225814	Battery Powered Decapper	13	1
W225820	Battery Powered Crimper	20	1
W225821	Battery Powered Decapper	20	1
W225815	Battery Powered Flip Cap Crimper	13	1
W225822	Battery Powered Flip Cap Crimper	20	1

## Automatic Crimping Tool Stand

Cat. No.	Description	Qty/Case
W225701	Crimping Tool Stand	1

## Stoppers and Seals

Cat. No.	Description	Qty/Case
<b>Stoppers</b>		
W224100-400	13 mm Ultra Pure straight plug stopper	1000
W224100-405	20 mm Ultra Pure straight plug stopper	1000
<b>Seals</b>		
224177-01	13 mm open top unlined aluminum seal, natural	1000
224182-01	13 mm center disc tear-out unlined aluminum seal, natural	1000
224192-01	13 mm complete tear-off unlined aluminum seal, natural	1000
224202	13 mm flip cap unlined aluminum seal, red	1000
W224207	13 mm flip cap tear-off unlined aluminum seal, red	1000
224178-01	20 mm open top unlined aluminum seal, natural	1000
224183-01	20 mm center disc tear-out unlined aluminum seal, natural	1000
224191	20 mm solid top unlined aluminum seal, natural	1000
224193-01	20 mm complete tear-off unlined aluminum seal, natural	1000
224203	20 mm flip cap unlined aluminum seal, red	1000
224208	20 mm flip cap tear-off unlined aluminum seal, red	1000



DWK Life Sciences GmbH  
 Otto-Schott-Str. 21  
 97877 Wertheim / Main  
 Germany

Phone: +49 (0) 6131 - 1445 4131  
 Fax: +49 (0) 6131 - 1445 4016  
 sales@DWK-LifeSciences.com  
 www.dwk.com